MANUAL DE USO • USER'S MANUAL • MANUALE D'USO

SIBEL SOUND DUO

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Technical Service: Tel. +34 93 433 54 50 e-mail: sat

COPYRIGHT

No part of this publication may be reproduced, transmitted, transcribed, stored in a back-up system or translated into any language or computer language in any form or by any means, electronic, mechanical, optical, chemical or manual without the express written consent from SIBEL S.A.U.

DISCLAIMER

SIBEL S.A.U. is responsible for the security, reliability and performance of this equipment only if:

- The place where the system is installed or used meets the requirements for electrical installations IEC and other applicable regulations.
- The repairs, revisions or modifications, in warranty or not, are made by technical staff **of SIBEL S.A.U.**
- The device is used by qualified staff in accordance with the recommendations stated in this User's Manual.

INDEX

1. SAFETY	4
1.1. INTENDED USE	4
1.2. DIRECTIONS FOR USE	4
1.3. USE LIMITATIONS. CONTRAINDICATIONS	5
1.4. ⚠ WARNINGS	5
1.5. DISPOSAL OF ELECTRICAL AND ELECTRONIC EQUIPMENT WASTE BY USERS IN PRIVATE HOUSEHOLDS IN THE EU	8
2. INSTALLATION INSTRUCTIONS	8
2.1. MAIN SPECIFICATIONS	9
2.2. LIST OF CONTENTS	10
2.3. CALIBRATION	11
3. INSTRUCTIONS FOR USE	11
3.1. LAYOUT OF CONTROLS AND CONNECTORS	11
3.2. USE OF THE EQUIPMENT	12
3.3. INSTALLATION OF THE BATTERIES	14
3.4. TEST CONFIGURATION	14
3.5. EQUIPMENT CONFIGURATION	15
3.6. TONAL TEST	16
3.7. DIAGNOSES	22
3.8. DATABASE	23
3.9. LINK WITH THE PC SOFTWARE	24
3.10. ACTIVATING OPTIONS	25
3.11. FIRMWARE UPDATE	29
3.12. TRANSPORT AND STORAGE	29
4. CLEANING AND MAINTENANCE	30
4.1. CLEANING	30
4.2. PREVENTATIVE MAINTENANCE	30
4.3. CORRECTIVE MAINTENANCE	30
5. TECHNICAL SPECIFICATIONS	31
6. APPLICABLE STANDARDS AND DIRECTIVES	33
7. SYMBOLS	34
ANNEX 1: ELECTROMAGNETIC COMPATIBILITY	35
ANNEX 2: COMPLIANCE WITH DATA PROTECTION REGULATIONS	37
ANNEX 3: CALCULATION OF DIAGNOSTICS	38

Thank you for choosing **SIBELSOUND DUO**

SIBELSOUND DUO has been developed by the RDI department of **SIBEL S.A.U.**

This product is manufactured in accordance with the European Sanitary Products Directive and with the best quality standards.

CE 0318 **PRODUCT IN COMPLIANCE WITH MEDICAL
DEVICE DIRECTIVE 93/42/EEC. CLASS IIa**

Revised
Date: 2019-11
Technical Director

Approved
Date: 2019-11
Sales Director

1. SAFETY

SIBELSOUND DUO is a 'screening' audiometer which allows users to perform audiometric explorations in a fast, simple, and reliable manner, for determining hearing thresholds both through air and bone conduction, and without masking.

The SIBELSOUND DUO audiometer has been designed to offer maximum security. Please, read all the instructions for use before using the equipment. Otherwise, this may result in user or patient injuries, or damage to the equipment and/or accessories.

For further technical information, or information about the product, please contact SIBEL, S.A.U. or your regular supplier.

1.1. INTENDED USE

The **SIBELSOUND DUO** audiometer is a desktop unit designed to diagnose and follow-up patients with hearing loss (hypoacusis) caused by noise, disease or trauma, by generating sound and vibration stimuli.

The following considerations must be noted.

- For the exclusive use of health professionals, supervised and/or trained by a doctor.
- For use in a health center or similar (unfit for outdoor use)
- Not intended for use in moving transport vehicles

1.2. DIRECTIONS FOR USE

Continuous exposure to high-level sounds, sound-intensity peaks, traumatic injuries, or some diseases may deteriorate the hearing function or cause hypoacusis. The **SIBELSOUND DUO** audiometer is an electromedical device designed to diagnose hearing loss by performing tonal tests, both through air and bone conduction (with and without masking, with narrow-band white noise). Optionally, automatic tests can be performed, based on the Hughson-Westlake protocol.

SIBELSOUND DUO is a desktop unit powered by 4 AA batteries or by an external power supply, and it includes air conduction headphones, a bone vibrator and a response button for the patient. The unit can connect to a PC and be used together with the W50 audiometry software in order to conduct real-time tests or to download data. Optionally, it can connect to a network via Ethernet, to allow for its integration into the Clinical Information System by exporting result files in standard formats (PDF, XML, CSV).

The audiometer is not designed for other conditions of use, nor for power sources other than those specified in this manual.

SIBELSOUND DUO was designed exclusively to be used by health professionals, under a doctor's supervision and training. Users are recommended to have been trained in audiometry techniques. Users should familiarize themselves with the operation of the equipment before using it with patients. All the information necessary for operation is available in this User Manual.

SIBELSOUND DUO is suitable for patients older than 4 years, with a physical and mental level that allows them to wear the headphones or the vibrator, and able to use the patient's response button or to give indications to the technician, using gestures, about their auditory perceptions.

1.3. USE LIMITATIONS. CONTRAINDICATIONS

The tests must be interpreted by a doctor, and only a doctor can indicate the treatment to be applied. The health professional must take into account the patient's symptoms before performing an audiometric test. The acceptability of a test is the responsibility of the health professional. The audiometer must not be used when the validity of the results might be compromised by external factors.

Audiometry tests require the patient's cooperation; the patient needs to press a warning button when they detect stimuli. The doctor needs to assess the patient's ability to undergo the audiometry test. Special attention must be paid to children, elderly people, and disable patients.

The EN ISO 8253-1 standard specifies certain requirements on the environmental noise conditions when performing audiometric tests. It is recommended to install the audiometer with a soundproof booth for the audiometric tests, or to adapt noise-canceling devices in the headphones, so that the validity of the results is not compromised. Please, refer to the booth user manual to connect the audiometer, or to the instructions for the noise canceling device.

1.4. WARNINGS

Use only the **SIBELSOUND DUO** with accessories provided by the manufacturer or dealer, or those that meet the specifications of this manual. The use of other accessories with the **SIBELSOUND DUO** can cause damage to the device, increase electromagnetic emissions or decrease immunity.

Handle accessories by their strongest parts (e.g. the connectors when manipulating the USB cable). **DO NOT** get them wet or expose them to very abrupt changes of temperature. **DO NOT** apply excessive stress to the accessories and the audiometer. In particular, avoid pulling or bending any part of the cables.

DO NOT disconnect the cable from the device by pulling the cable. You can damage the device or accessories reducing the product's safety. Hold always the cable by the connector.

Deteriorated accessories (power supply or USB cable and accessories with broken cases, cables or connectors) should not be used since there is risk of cramps.

Please contact with [SIBEL S.A.U.](#) or your dealer to get new accessories.

In case of receiving the product with a deteriorated packaging, **DO NOT** use the device and please contact your courier agency, dealer or [SIBEL S.A.U.](#) After Sales Service.

DO NOT use the **SIBELSOUND DUO** out of the specified environmental conditions of use and storage.

DO NOT use the **SIBELSOUND DUO** in moving transport vehicles.

DO NOT use transducers (headphones or bone vibrator) that have not been calibrated with the audiometer. The label of each accessory indicates the serial number of the audiometer with which it has been calibrated.

In case of device over-heating during use, immediately turn off the device and please contact [SIBEL's After Sales Service](#) or your dealer.

The equipment should be placed in a safe position to avoid falling that may result in equipment damage or harm to the patient and/or user.

Keep your device protected from shock and vibration. During transportation, place all the items inside the carrying bag. It provides enough protection against small accidental impact. Always transport the device and accessories inside the carrying bag. It provides protection against small accidental impacts.

Remove the batteries if you do not use the audiometer for a long period of time.

ELECTRICAL RISK

The **SIBELSOUND DUO** is not intended to be used with other energy sources that are not covered in this manual.

DO NOT try to open the device. In case of malfunction or unexpected operation, please contact [SIBEL's After Sales Service](#) or your dealer.

DO NOT use multiple mains sockets to not degrade the electrical safety.

Following standard EN60601-1, the unit is classified as continuous operation mode.

Essential operation includes generating acoustic or vibratory signals (complying with the specified accuracy and/or ensuring that it is not possible to make an incorrect diagnosis of the patient), presenting hearing levels, and storing the complete data into the memory. Communication with a computer is not considered essential operation, provided that the storage of the complete data is guaranteed.

RISK OF ELECTRICAL SHOCK

No parts are allowed for temporary immersion. MAY CAUSE ELECTRIC DISCHARGE

Contact of liquids with the internal parts of the device and the connectors must always be avoided.

SIBELSOUND DUO is designed to be used in conjunction with a software application such as the [W50](#) Audiometry Software, which runs on an external device. This device must comply with the current standards of electrical safety and electromagnetic compatibility.

RISK OF EXPLOSION

The system is not designed to work in an explosive environment or in the presence of flammable anesthetics or gases of any kind. MAY CAUSE EXPLOSION.

RISK OF CONTAMINATION

Perform cleaning and maintenance of device and accessories with the specified frequency and following the instructions in this manual in order to keep product's safety.

RISK OF INTERFERENCE

This is an electronic medical device; therefore requires special precautions regarding electromagnetic compatibility (EMC): it must be installed and put in service according to the information attached in [Annex 4. ELECTROMAGNETIC COMPATIBILITY](#).

DO NOT use the system in the presence of radio equipment, such as mobile phones, transmitters and similar equipment generating radio frequency emissions. Follow the recommendations regarding the separation distance specified in the manufacturer's declaration on electromagnetic compatibility in this manual.

DO NOT use the **SIBELSOUND DUO** adjacent or stacked with other equipment. It could result in improper operation.

DO NOT use the system in and MRI environment.

1.5. DISPOSAL OF ELECTRICAL AND ELECTRONIC EQUIPMENT WASTE BY USERS IN PRIVATE HOUSEHOLDS IN THE EU

  **DO NOT** dispose **SIBELSOUND DUO** audiometer, accessories and / or batteries with household waste. Deliver it in a designated collection point for recycling in accordance with the legal requirements of your country.

The unit uses a lithium battery, and may use NiMh or Alkaline batteries.

To obtain further information on the correct disposal of the product and its accessories, please contact **SIBEL S.A.U.** Customer Service or your distributor.

2. INSTALLATION INSTRUCTIONS

Before using the unit for the first time, select the language as described in section 3.5 EQUIPMENT CONFIGURATION.

Connect all the necessary accessories to the back of the unit, as described in section 3.1, ensuring that the cables do not interfere with the work.

2.1. MAIN SPECIFICATIONS

SIBELSOUND DUO is available in two models, with the following specifications:

FUNCTIONALITIES	A	AOM
Two Channels		
Pure Tone Air Conduction		
Pure Tone Bone Conduction		
Masking with narrow band noise		
Intermediate Frequencies (125, 750, 1500 Hz)		
Internal Data Base (> 1000 tests)		
Internal Data Base (> 1000 tests)		
Frequency modulation (Warble)		
Diagnosis (Coulcil, IMSS Mexico, ELI, SAL, Klockhoff)		
Calculation of hearing loss		
Audiometry Software W50 (Demo)		

Intercommunicator		
Firmware update		

CONNECTIVITY		
USB 2.0 to PC		
Ethernet		
OPTIONS		
Hugshon-Westlake automatic test		
Data export (XML, PDF or CSV)		
Audiometry Software W50 (Licence)		

Included

Optional

Not included

This audiometer has been manufactured following strict quality controls. However, accidents may occur during transportation or storage. For this reason, it is advisable to check the state of the unit and its accessories before installation. Do not discard the packaging until correct operation of the unit has been completely verified.

2.2. LIST OF CONTENTS

COD.	QUANTITY		DESCRIPTION	A	AOM
AUDIOMETER					
---	1		SIBELSOUND DUO SN: 208 - _____		
---	1		USER'S MANUAL SIBELSOUND DUO (DOC.520-800-MU__REV: __)		

STANDARD ACCESSORIES					
03658	1		USB CABLE TYPE A-B 2.0, PC CONNECTION		
03024	1		DEMO MODE W50 AUDIOMETRY SOFTWARE		
01807	1		PATIENT RESPONSE SWITCH		
09460	1		AIR TRANSDUCER SET (DD45)		
02258	1		SET OF HEADPHONES, BONE CONDUCTION B71W		
09489	1		CARRYING BAG		
09459	1		SIBELSOUND DUO EXTERNAL POWER SUPPLY		

OPTIONAL ACCESSORIES					
09134	1	-	AUTOMATIC AUDIOMETRY OPTION		
09505	1		ETHERNET EXPORT OPTION • ETHERNET CABLE • FIRMWARE OPTION		
03673	4		ALCALINE BATTERY 1.5V LR6 A		
02303	1		W50 SOFTWARE LICENSE		
02260	1		AIR CONDUCTION SET (TDH39)		
03078	1		SOUND SUPPRESSORS FOR AIR CONDUCTION (AUDIOCUPS)		

Included

Optional

Not ttt

2.3. CALIBRATION

It is recommended to calibrate the transducers and to perform a general check of their safety systems, settings, functions, etc. with a frequency not exceeding twelve months (EN ISO 8253-1), or when malfunction of the unit is suspected. The unit must be calibrated following the manufacturer's indications, available for the Technical Departments or personnel authorized by [SIBEL S.A.U.](#)

3. INSTRUCTIONS FOR USE

3.1. LAYOUT OF CONTROLS AND CONNECTORS

- A** ON / OFF Button
- B** Intensity adjustment
- C** Frequency adjustment
- D** Signal button
- E** Touch screen
- F** Patient response connector
- G** Bone transducer connector
- H** Right headphone connector
- I** Left headphone connector
- J** Ethernet connector
- K** USB connector
- L** External power connector


3.2. USE OF THE EQUIPMENT

To switch on the unit, press the ON/OFF button (A) on the lower part of the screen. During start-up, the unit performs a self-check process for all connected accessories and for calibration. If a problem is detected, a warning screen appears with the non-connected accessories and/or the last calibration date.

MAIN SCREEN

- a** Date and time setting and indication
- b** Identifier for the current patient and access to editing of the patient's details
- c** Ethernet connection indicator
- d** Screen brightness setting
- e** Battery / power supply indicator and selection of the type of battery
- f** Continual / Pulsating tone
- g** Pure / Warble Tone (FM)
- h** Masking level selection and indicator
- i** Selection of the ear to be checked, Right (R) / Left (L)
- j** Intensity indicator
- k** Frequency indicator

Main screen items can be customized so that the unit always starts up with the same configuration (please, refer to section 3.5 EQUIPMENT CONFIGURATION).

The main screen allows users to directly and immediately perform a test. When pressing the *Sound Application* button (D), the audiometer will emit the selected tone, in the selected mode and ear. On the center of the screen, the symbol  will appear to indicate that the signal is being applied to the patient. To change the signal application mode to Inverted Mode, press the C and D buttons simultaneously. Press again to deactivate Inverted Mode.

Press the **C** button to activate the intercom, and keep pressing it while talking to the patient. On the center of the screen, the icon  will indicate that the intercom is in use.

When pressing  **Mask**, the masking value will flash; turn the intensity setting (**B**) to change it.

Press  **A** to change the test mode, between air (A), air with masking (AM), bone (O), or bone with masking (OM). The A model only offers the air mode (A). If the Automatic Audiometry option is available, this test mode can also be selected (Auto HW).

Press  to display the data for the current test.

Press  to create a new patient.

Press  to save the current test (the icon changes to  indicate that the test has already been stored).

Press  to access the Main Menu, with the following options:

-  Test configuration (refer to section 3.4)
-  Equipment configuration (refer to section 3.5)
-  Database (refer to section 3.8)
-  Maintenance (refer to section 4)
-  Back
-  Main Screen

3.3. INSTALLATION OF THE BATTERIES

 **DO NOT use battery types other than the ones specified in this manual. DO NOT use deteriorated batteries since there is risk of skin irritation and sensitization.**

 **Respect the polarity indicated in the battery compartment. The device will not turn on otherwise.**

To install or change the batteries in the **SIBELSOUND DUO**, follow these steps:

- Switch the unit off. If you have an external power supply, make sure it is not connected to the unit.
- Open the battery compartment loosening the screw.
- Remove the used batteries and/or insert new ones following polarity shown.
- Close the battery compartment re-tightening the screw.
- Switch the unit on again (or connect the power supply).

SIBELSOUND DUO has a battery status indicator (e). If the indicator light is yellow , this means that there is less than 25% battery charge; if it is red , the battery charge is below 10%. Replace batteries as soon as possible.

To optimize battery life, it is recommended to set brightness to 50%

When the unit is not used for longer than 3 minutes, the screen will switch off to save energy. Touch the screen or any button to reactivate the screen.

3.4. TEST CONFIGURATION

This menu allows users to customize audiometric tests in all modes, depending on their needs.

Press  to select frequencies to be explored during tonal tests.

Press  to select the different diagnoses available (see section 3.7).

Press  to personalize the main screen default values, which will apply every time the unit starts up (mode, tone type, starting ear, intensity and initial frequencies, etc.).

Press  to select the threshold saving method (manual, automatic as per the latest level applied, or automatic as per the patient's response).

Press  to activate automatic masking. When activated, the Main Screen displays the  icon, and the masking level changes automatically when the applied intensity is modified, so that the difference remains constant. When deactivated, the Main Screen displays the  icon, and the masking level needs to be modified manually and independently from the selected signal level.

 indicates that the HL correction is being applied, as per the ANSI standard. Press to change the type of correction.

 indicates that the HL correction is being applied, as per the ISO standard.

 indicates that the Inverted Tone is disabled. Press to activate;  will display on the screen.

Press  to access the automatic Hughson-Westlake audiometry test configuration, optional in all models (refer to section 3.10).

3.5. EQUIPMENT CONFIGURATION

This menu allows users to customize the configuration of the unit according to their needs.

Press  to select the language of the unit.

Press  to set date and time. This configuration can also be accessed pressing the date/time indicator  on the top bar of any screen.

Press  to set the configuration template stored in the unit.

Press  to store the current equipment configuration as Template Configuration, so that it can be restored at any time.

Press  to protect access to the unit using a PIN code.

Press  to access the Ethernet connection and data export settings, optional in all models (refer to section 3.10).

3.6. TONAL TEST

GENERAL CONSIDERATIONS

Tonal audiometry is the basic test for the diagnosis of hearing loss; it consists of determining hearing threshold levels via air and/or bone conduction, with or without masking.

There are different procedures for conducting the test; the health professional must determine which procedure is most appropriate for each patient. It is recommended to explore, at least, 4 frequencies through air conduction (500, 1000, 2000 and 3000 or 4000 Hz) for both ears.

Tests should be carried out by a technician trained on how to perform audiometric tests, since they need to decide which ear is to be tested first (usually, the one with better hearing), whether or not masking is required, if the patient is producing coherent responses (**SIBEL SOUND DUO** warns in case of discrepancies) and whether or not it is necessary to interrupt the test.

Preferably, tests should last less than 20 min, to prevent patient fatigue.

PREPARING AND INSTRUCTING THE PATIENT

Before performing the test, perform an otoscopic examination of the patient in order to identify possible obstructions in the external auditory canal which would need extracting.

It is recommended that the patient arrives at the test site at least 5 minutes in advance, to avoid possible exposure to high-noise levels, or previous physical activities which might alter the test results.

Before starting, the technician should give the patient instructions about the test procedure, indicating, amongst others, the following:

- Which ear will be examined first.
- The general tone presentation sequence.
- How to indicate a response using the patient's response button (or gestures).
- When the patient is expected to provide a response (when they hear the tone, not when they think they hear it).

Patients need to take off their glasses and/or accessories, and tie their hair back, whenever these may prevent the correct placement of the headphones or the bone vibrator. When placing the air conduction headphones, make sure that the **RED** headphone is on the patient's right ear, and the **BLUE** one on the left ear.

PATIENT'S DETAILS

Before starting the test, press  on the main screen to register the patient in the audiometer. The patient data editing screen will be displayed; enter their identification data (ID number, clinical history number, etc.), name and surname. If the ELI diagnosis is selected, the patient's age and gender also need to be entered.

Optionally, the automatic Identifier mode can also be activated. The audiometer will automatically assign a sequential number to the new patient (preceded by the # symbol in the identifier), so the patient's name and surname will not need to be entered.

When this mode is activated, pressing  on the main screen creates a new patient with the appropriate identifier number, and the patient's data editing screen does not come up.

A previously-saved patient can also be selected from the database by pressing  on the patient's data screen.

At any time during the test, patient data can be accessed by pressing the patient's identifier box (b) on the top bar of the screen. Note: When creating a new patient, if there is an ongoing test with unsaved values, a warning message displays, so that these can be saved before starting a new test.

SETTING A HEARING THRESHOLD

To set the hearing threshold at a certain frequency, follow these steps:

- Select the frequency that you want to examine by turning the Frequency Setting (C). Normally, this starts at the 1000 Hz frequency.
- Select the signal level you want to apply by turning the Intensity Setting (B). You can start with a level which is clearly audible for healthy people, for example, 40 dB.
- Apply a continuous tone for 1 or 2 seconds, pressing the Signal Application button (D).
- When the patient presses their Response Button, the main screen background will turn green.
- If the patient responds before the signal is applied, or more than 2 seconds after it is applied, the screen background will turn red, indicating that the patient's response is not coherent. Apply the tone again until you obtain a coherent response, or give the patient instructions again.
- If the patient gives a correct response, reduce the intensity by 10 dB and reapply the tone.
- Continue decreasing intensity in 10 dB steps and applying tones until there is no response from the patient.
- If the patient does not respond, increase intensity in 5 dB steps, until they respond.
- Every time the patient produces a response, decrease the intensity by 10 dB, and every time there is no response, increase it by 5 dB.
- The time interval between tones must be varied randomly, so that the patient does not predict when the tone is going to be applied. The interval between two tones must always be longer than their duration.

- Once you have obtained 2 responses for 3 attempts at the same level (or 3 responses for 5 attempts), this will be the hearing threshold for the selected frequency.
- Press the Intensity Setting (B) to store the threshold value. The intensity and frequency values will be displayed in green for 1 second.

The audiometer can be configured so that it stores thresholds automatically, without having to press the Intensity Setting.

To activate this option, access the TEST CONFIGURATION  menu and press . A menu will be displayed to select between two automatic saving methods:

- Last level applied: the last intensity level for which a signal has been applied will be saved as the threshold.
- Patient's response: the last level applied for which the patient has given a correct response using the response button will be saved as the threshold.

To perform an air conduction tonal test, follow these steps:

- On the main screen, select the air conduction mode (A).
- Select the first ear to be examined (R/L), the tone presentation mode (continuous  or pulsating ) and the type of tone (Pure  or Warble .
- Select the 1000 Hz frequency, set the hearing threshold, and store it following the steps described above.
- Then, select the frequency immediately above, and set the appropriate hearing threshold.
- Continue examining the different frequencies in ascending order, until reaching the highest one. Then, select the frequency immediately below 1000 Hz, and set the hearing thresholds in descending order (for example: 1000, 1500, 2000, 3000, 4000, 6000, 8000, 750, 500, 250, 125).
- It is not necessary to examine all the frequencies. The audiometer can be configured so that it only shows those you intend to examine. This is done through the test configuration menu (refer to section 3.4).
- Once the first ear has been examined, press R (or L) to change to the other one, and examine all the frequencies for the second ear.
- At any time during the test, you can view the stored thresholds by pressing  on the main screen (refer to section TEST DATA).
- Once all thresholds have been set (or at any time during the test), press  on the main screen (or on the Test Data screen) to save it. If the Data Export Option is available, and Automatic Export is on, the files will be automatically sent to your network.

BONE CONDUCTION TONE TEST

The bone vibrator must be handled with care because small hits may damage it.

In this type of test, the tone is applied through a bone vibrator. It is very important to position the vibrator correctly so that the results are reliable. The vibrator is placed on the mastoid bone; the tone applied will be a few decibels above the hearing threshold, instruct the patient to move it along the mastoid bone until they find the point where they hear it the most. Make sure that the vibrator is perfectly attached to the mastoid bone, and that it does not touch the auricle, to prevent cartilaginous conduction. The bone-conduction tonal test is performed exactly in the same way as the air-conduction one, except for the following:

- The bone vibrator must be placed on the patient's mastoid bone, on the first ear to be examined (do not use the headphones).
- On the main screen, select the bone conduction mode ().
- The bone-conduction examination does not allow 125 and 8000 Hz frequencies.
- Before examining the second ear, select it on the main screen and change the position of the patient's bone vibrator.

AIR-CONDUCTION TONAL TEST WITH MASKING

This type of test is performed to prevent the patient from hearing the test tone through the ear opposite to the one being examined; narrow-band noise is applied on that ear. It is recommended to follow the steps below to set the hearing threshold with masking:

- On the main screen, select the air conduction with masking mode ().
- Select the ear to be examined and set the intensity to a level equal to the hearing threshold without masking.
- Press  on the main screen. The dB value indicated in this icon will flash.
- Turn the Intensity Setting () to set the initial masking level to the hearing threshold value of the masked ear.
- Press  again or press the Intensity Setting () to save the masking level (the value no longer flashes)
- Press the Signal Application button () to apply the tone. The tone will be applied on the ear to be examined, and the masking noise to the other one.
- Increase the masking level until the tone is inaudible, or until it exceeds the tone level applied on the ear to be examined.
- If the tone is still audible when the masking level is at the tone level, this will be the hearing threshold. Press the Intensity Setting () to store it. The intensity and frequency values will be displayed in green for 1 second.
- If the tone is masked, increase its level until it is audible again.
- Increase the masking level by 5 dB. If the tone is inaudible, increase its level until it is audible again.
- Repeat this process until the tone remains audible, even if the masking level has been increased more than 10 dB. This will be the threshold of hearing. Press the Intensity Setting () to store it.

- Repeat this process until all the desired frequencies have been examined, following the same order as in the tonal test without masking.

If you want the difference between the masking level and the applied tone to always be the same, activate Automatic Masking mode through the Test Configuration menu (refer to section 3.4).

- On the main screen, adjust the desired tone level.
- Then, adjust the masking level by pressing  and turning the Intensity Setting (B), until the desired difference is obtained.
- Press  again or press the Intensity Setting (B) to save the masking level.
- From now on, when the tone level is modified, the masking level will change automatically in order to maintain a constant difference.

BONE-CONDUCTION TONAL TEST WITH MASKING

The bone-conduction tonal test with masking is performed exactly in the same way as the air-conduction test with masking, except for the following:

- Place the headphones on both of the patient's ears, and the bone vibrator on the mastoid bone of the first ear to be examined (check the recommendations for placing the bone vibrator correctly in section [BONE-CONDUCTION TONAL TEST](#)).
- On the main screen, select the bone conduction mode with masking (.
- When increasing the masking level for the first time, do it until the tone is inaudible, or until it exceeds by 40 dB the tone level applied to the ear to be examined.
- If the tone is still audible when the masking level is 40 dB over the tone level, this will be the hearing threshold.
- Continue with the test in the same way as with the bone-conduction test with masking, until the tone remains audible even if the masking level has been increased more than 10 dB. This will be the hearing threshold.
- Examine the other desired frequencies (125 and 8000 Hz frequencies not allowed).
- Before examining the second ear, select it on the main screen and only change the position of the patient's bone vibrator, keeping the headphones on both ears.

TEST DATA

Press  on the main screen to display a data table for the current test similar to the one below:

The rows indicate the ear and the type of test:

- RA: right ear, air conduction
- LA: left ear, air conduction
- RO: right ear, bone conduction
- LO: left ear, bone conduction

The columns show the examined frequencies with the threshold values saved (in dB). Turn the Frequency Setting (D) to move the table columns to the left and right. To modify any of the table values, press the value to display the main window, with the corresponding intensity and frequency, in order to redefine such threshold.

Note: if the data displayed refer to an automatic test (Hughson-Westlake), this will be indicated in the top left corner of the table with .

27/06/19 11:15

ID: AB0123456789Y

	500	750	1000	1500	2000	3000	4000
RA (M)	5 (20)	5 (20)	10 (25)	15 (30)	35 (50)	30 (45)	15 (30)
LA	0 (15)	5 (20)	0 (15)	5 (20)	0 (15)	0 (15)	5 (20)
RO	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)
LO	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)	--- (---)

Press  to view the masking values, displayed in brackets below each threshold values. If the test is performed without masking, this button will be grayed out.

Press  to display the results of the selected diagnoses (see section 3.7)

Press  to save the test.

3.7. DIAGNOSES

The **SIBEL SOUND DUO** audiometer calculates the following diagnoses:

- 'MINISTRY' DIAGNOSIS, IN ACCORDANCE WITH THE MINISTRY OF LABOR AND SOCIAL AFFAIRS, BOE (OFFICIAL STATE BULLETIN) NO. 22 (JAN 26 2000).
- COUNCIL OF PHYSICAL THERAPY DIAGNOSIS
- DIAGNOSIS IN ACCORDANCE WITH THE MEXICAN SOCIAL SECURITY INSTITUTE.
- ELI INDEX (Early Loss Index).
- SAL INDEX (Speech Average Loss Index).
- KLOCKHOOF DIAGNOSIS (Modified by Lavoro Clinic in Milan).

From the Test Data screen, the user can

press  to access the diagnoses display. The following screen will appear:

For every diagnosis, there is a table with values for Hearing Loss, Average dBs and Bilateral Loss in each ear. Then, the diagnostic classification is displayed for each ear. For the ELI and Klockhoff indexes, only the diagnostic classification for each ear is displayed (without hearing loss values).

Turn the Intensity Setting (B) to scroll the screen up or down, and view all the diagnoses. Only diagnoses selected in the menu TEST CONFIGURATION  -> DIAGNOSES  will be displayed.

Refer to [Annex 3](#) to obtain information about calculation methods for each diagnosis.

27/06/19 11:15

ID: AB0123456789Y

COUNCIL OF PHYSICAL THERAPY	R.E.	L.E.
Hearing Loss (%)	14,1	2,7
Average dBs	26,2	13,8
Bilateral Loss (%)	4,1	
Right Ear	Slight Hypoacusia	
Left Ear	Normal	

3.8. DATABASE

The SIBELSOUND DUO audiometer has a database for up to 1024 tests. Every time a test is saved by pressing , it is automatically stored in the unit's internal database. Stored records can be viewed in

the database by pressing  from the Main Menu. The following screen appears:

The records table consists of 5 columns with the following information:

- Record selection marks
- ID: patient's ID
- Name: patient's full name
- Date-Time: test start date and time
- Exp: whether or not the record has been exported

Turn the Intensity Setting (B) to vertically scroll the rows cursor (dark blue) through the database records. Press each column heading to sort the records in ascending or descending order, depending on the field. From the database screen, the following actions can be performed:

27/06/2019 11:15

ID:AB0123456789Y

☐	ID	Name	Date-Time ▼	Exp
☐	42386594N	PÉREZ CAPDEVILA, ...	27/06/19-09:59	-
☒	27456381K	DEL HORNO GARCÍ...	22/06/19-11:03	-
☒	537621984	MORENO, FRANCIS..	22/06/19-12:31	-
☒	52859365Q	DÍAZ, JOSÉ	20/06/19-15:25	-
☐	85472843J	ORTIZ, MANUEL	15/06/19-09:03	-
☐	27456381K	DEL HORNO GARCÍ...	17/06/19-09:45	-
☐	42386594N	PÉREZ CAPDEVILA, ...	14/06/19-10:41	-

- Press  to search the database by text or date.
- If the Ethernet / Data Export option is available, one or several records can be exported to an external network (refer to section 3.10). First, select the records in question by pressing the box for each row , or placing the rows cursor on the desired record and pressing the Intensity Setting button (B). It is also possible to select all the database records by pressing the box located on the heading of the table's first column.
- Press  to export the selected records. A message will specify the number of records to be exported and, once exported correctly, this will be indicated with the text 'OK' in column 'Exp'.
- Press  to delete selected records.

- To delete all the database contents, press . Two warning messages appear; they need to be confirmed before proceeding, since this action will COMPLETELY DELETE ALL THE INFORMATION in the database.
- To start a new test with a patient already stored in the database, choose the patient by selecting the right-hand-side box, and then press . The Home Screen appears, with the selected patient record ID.

The  button only appears active when there is just one record selected. If there are no selected records, or if there is more than one, this button will be grayed out.

To review the contents of each record, place the rows cursor on the desired record, and press . The Test Data table appears, where the user can view the stored thresholds, and the diagnoses results.

3.9. LINK WITH THE PC SOFTWARE

The **SIBELSOUND DUO** audiometer can be used together with the [W50 Audiometry Software](#), to perform real-time audiometric tests, viewing the data directly on the PC screen. To do this, connect a personal computer using a USB cable to the appropriate connector on the back panel (K).

Real-time tests must be performed from the audiometer's main screen. The software manages processes such as patient's registration, threshold storage and test storage. Consequently, the following main screen features are disabled when a real-time test is initiated:

- Test data
- New patient
- Save test

Real-time tests are not stored in the database, and no diagnoses are calculated.

From the [W50](#) software, the user can also download the contents of their **SIBELSOUND DUO** database, or the calibration record in CSV format.

For information about the procedure for real-time tests, and about downloading information to a PC, please refer to the [W50](#) software user manual.

3.10. ACTIVATING OPTIONS

Both **SIBELSOUND DUO** models have two additional features which can be optionally obtained in order to expand the unit's capabilities:

- Ethernet / Data Export option
- Automatic Audiometry (Hughson-Westlake) option

To enable these options, enter an activation code in the unit. To do this, access the Maintenance menu

() and press () . A screen appears with 4 fields for the code, the unit's serial number and a numeric keypad.

Enter the code and press () . If the code is correct, a message appears confirming that the option is active. Repeat this process to enable the other option.

ETHERNET / DATA EXPORT OPTION

With the ETHERNET / DATA EXPORT option, test results can be transferred to a network in three different formats (PDF, XML or CSV), for further review and analysis.

To configure the SIBELSOUND DUO Ethernet connection and data export, access the EQUIPMENT CONFIGURATION menu

() and press () . The Ethernet configuration should ideally be carried out by an IT technician who is familiar with the network where the audiometer is being connected.

Press () to configure the IPv4 network. A screen appears where the user can introduce manually the network data for IP addresses, mask, gateway, primary DNS and secondary DNS. Optionally, the Auto IP mode can be selected, so that the network assigns these addresses for the unit.

Press () to set the User, Password and Domain for network access.

Press () to select the format type for the files to be exported (PDF, XML or CSV) and the destination folders in your network. The Auto Export mode can also be activated, so that every time a test is saved, the current test files are automatically exported to your network. When records are exported from the database, the selected files will be generated in this screen, and saved to the same folders.

If your computer runs on Windows 10, compatibility with the Samba (SMB) needs to be enabled. To do this, access Control Panel / Programs and Features / Enable or Disable Windows Features, or check with your IT technician.

If there is an error during the export process, a warning message appears indicating the type of error.

For more information about export files format, please contact [SIBEL S.A.U.](#) Customer Service.

AUTOMATIC AUDIOMETRY (HUGHSON-WESTLAKE) OPTION


The AUTOMATIC AUDIOMETRY (HUGHSON-WESTLAKE) option allows users to perform an air-conduction tonal test (without masking) in a fast, simple manner, with minimum technician intervention.

To personalize various parameters for the automatic test, access the TEST CONFIGURATION menu  and press . A screen appears, where the following automatic test items can be customized:

- Threshold method (2/3 or 3/5)
- Tone type (Continuous or Pulsating)
- Starting level (between 0 and 70 dB)
- Random Off time (between 0 and 2 seconds)
- Tone duration (between 0.5 and 2 seconds)
- Frequencies to examine (125, 250, 500, 750, 1500, 2000, 3000, 4000, 6000, 8000 Hz). The 1000 Hz frequency is always examined, and does not appear on this screen.

To perform an automatic test, select  mode on the main screen. The information displayed depends upon:

- l** Tone type (depending on the test configuration selections)
- m** Test start (Play)
- n** Test stop (Stop)
- o** Familiarization
- p** Selecting the first ear to examine, Right (R) / Left (L)

The automatic test is always performed with the autosave by patient's response mode enabled. The procedure is as follows:

- Press  to create a new patient (see section [PATIENT'S DETAILS](#) in section 3.6).
- Press  or  to select the first ear to examine.
- Press  to initiate the familiarization process and instruct the patient on the test procedure. This starts with 1000 Hz tones at different intensities, which the patient needs to respond to by pressing the Patient's Response Button. This process identifies a level of intensity close to the threshold, which is used to start the automatic test.
- After the familiarization process, press  to initiate the automatic test. The unit will issue different tones at different frequencies, following a process similar to the one in the air-conduction tonal test. It will start by issuing a 1000 Hz tone with an intensity 10 dB above the one identified during the familiarization stage, or with the intensity set as Starting Level on the automatic test configuration screen.

- The duration and type of tones are those defined on the automatic test configuration screen. The tone separation intervals are random, with a minimum duration of 2 seconds and a maximum one of 2 seconds plus the random Off time set in the configuration process.
- After determining the frequency threshold, the Intensity (j) and Frequency (k) values turn green for a second, and the threshold is stored.
- The audiometer automatically changes to the following frequency selected in the automatic test configuration screen, following the same order as for the tonal tests. The first tone level applied at each frequency will be the same as for the 1000 Hz frequency.
- After the high frequencies have been examined, the audiometer examines again the 1000 Hz frequency to check the consistency of the results. If the difference between the first and the second 1000 Hz exam is equal or more than 10 dB, the message 'Inconsistent results' appears, indicating the difference value, and the test stops. If the difference is below 10 dB, the test continues, examining the lower frequencies.
- Once all frequencies for the first ear have been examined, the audiometer automatically changes to the other ear and examines it for all the frequencies. For the second ear, the 1000 Hz examination is only performed once.
- After the process, the message 'Test completed' appears.
- Press  to save the test. If the Data Export Option is available, and Automatic Export is on, the files will be automatically sent to your network.
- Press  to display the saved thresholds and diagnoses.

At any point, press  to pause the test. Press  again to resume the test. If the threshold for the current frequency has not yet been set when pausing the test, the test restarts at the same frequency, but with the Familiarization level or the Starting Level as set during Configuration. After pausing, press  to display the test data.

Press  at any time to stop the test. This concludes the current test. The user can save the test and view the results but, if they want to resume it, a new test needs to be initiated.

If you want to speak to the patient during the test, press the Frequency Setting (D) to activate the intercom. This pauses the test as when  is pressed.

3.11. FIRMWARE UPDATE

The **SIBELSOUND DUO** audiometer includes a feature to update the unit's program, and include all the improvements that **SIBEL S.A.U.** constantly develops for its products. Follow the steps below to update the firmware:

- Ensure that the unit is on and correctly connected to the computer via the USB cable.
- Press the Intensity Setting (**B**) 7 times and the Signal button (**D**) once to activate the firmware update mode.
- A start-up screen appears, while the system waits for the update file.
- From **W50 Software**, run the operation to send the update file (refer to the **W50 Software** user manual).
- After checking that the file sent is correct, the unit updates the firmware automatically.
- Finally, the unit indicates if the firmware has been successfully updated.
- If an error occurs during the update process, the unit informs the user, and the program is not modified.
- Once the update process has been completed, press the screen to restart the unit.

Note: Do not cancel the firmware update process once started.

3.12. TRANSPORT AND STORAGE

The two **SIBELSOUND DUO** models include a transport bag to store the unit and its accessories, for easy transportation. This bag can also be used to store the unit when not in use. It is recommended to remove the batteries from the battery compartment if the unit is not going to be used for over a week.

4. CLEANING AND MAINTENANCE

4.1. CLEANING

Do not use abrasive substances or solvents. The device and accessories can be damaged.

Clean the unit and its accessories with a wet, soft cloth, with water and a mild soap, or with isopropyl alcohol. Then, dry all the parts thoroughly. Ensure that no liquid or foreign body enters the unit, connectors or connections. Exercise caution not to scratch the screen.

The air-conduction headphone pads can be removed for thorough cleaning. Dry them with a cloth before replacing.

4.2. PREVENTATIVE MAINTENANCE

Preventative maintenance consists of all those actions that keep the equipment in good working order. It is recommended to do the following every week:

- Verify that all connections, accessories and other external elements are in perfect condition, without cracks or external damage.
- Check that all accessories are properly connected to the unit's back panel.
- Check the correct calibration of the unit. If a calibration warning appears when starting the unit, contact your supplier or [SIBEL S.A.U. Technical Department](#).

Some unit basic checks can be performed accessing the MAINTENANCE () -> UNIT TEST () menu and following the on-screen instructions.

If some anomaly or unusual unit behavior is detected, contact your supplier or [SIBEL S.A.U. Technical Department](#).

4.3. CORRECTIVE MAINTENANCE

The corrective maintenance consists of repairing and restoring the use of the equipment, when poor operation or improper use may have affected its functionality.

If some fault is detected, which prevents the correct operation of the unit, contact your technical support service or [SIBEL, S.A.U. Technical Department](#), indicating the type of anomaly encountered. The technical documentation necessary to perform the unit's corrective maintenance is available for the technical support services or for technical staff authorized by [SIBEL S.A.U.](#)

5. TECHNICAL SPECIFICATIONS

Healthcare Product Classification	Class IIa
Electrical Safety Classification (EN 60601-1)	Class II / Internally powered Type B applicable parts
Audiometer Type Classification (EN 60645-1)	Model A: Type 4. AOM model: Type 3.
Graphic Touch Screen	95 x 54 (mm.) 480 x 272 (pix.)
ON-OFF Button	Capacitive
Weight (g.)	600 approx. (without batteries and accessories)
Dimensions (mm.)	236 x 197 x 77
Connectivity	USB 2.0 (Communication with W50 Software) Ethernet (Data export via the network)
Power Supply	4 AA 1.5 V batteries, or external power supply (model: ACM06US12 o UES06Wx-120050SPA): input 100-240 V AC, output: 12 V DC
Battery Life	7 hours approximately, depending on use
Operating Conditions	Temperature 0 to 40 °C Humidity: 15% to 90% (without condensation) Pressure: 850 to 1060 hPa.
Transport and Storage	Temperature -20 to 70 °C Humidity: < 93 % (without condensation) Pressure: 850 to 1060 hPa.
Useful Life	7 years
Database	> 1,000 tests
Test Types	Tonal audiometry. Manual or Automatic. 2 channels. Automatic Hughson-Westlake audiometry.
Transducers	Air conduction: RadioEar DD45 or Telephonics TDH39 headphones. Reference thresholds: ISO 389-1 and EN 60318-1 - Bone conduction: RadioEar B71W vibrator. Reference thresholds: ISO 389-3. - Nominal impedance for headphones and vibrator: 10 Ohm.

Signal Sources	- Pure tones (see frequencies in the Table below). - Narrow-band noise (NBN). Band as per EN 60645-1. - Frequency Modulation (FM): Modulation: sinusoidal – 5Hz. Modulation percentage: 10%.
Type of Signal	Continuous and pulsating.
Signal Increase	5dB
Accuracy	Sound pressure: +/- 3 dB Vibratory force: +/- 5 dB
Linearity	+/- 1 dB
Frequency Accuracy	1%

Air Conduction Signal Levels

Hz	125	250	500	750	1000	1500	2000	3000	4000	6000	8000
Max. dB HL	70	85	100	100	100	100	100	100	100	100	85
Min. dB HL	-10	-10	-10	-10	-10	-10	-10	-10	-10	-10	-10

Masking Signal Levels (NBN)

Hz	125	250	500	750	1000	1500	2000	3000	4000	6000	8000
Max. dB HL	70	85	100	100	100	100	100	100	100	95	85
Min. dB HL	-10	-10	-10	-10	-10	-10	-10	-10	-10	-10	-10

Air Conduction Signal Levels

Hz	250	500	750	1000	1500	2000	3000	4000	6000
Max. dB HL	45	60	60	65	65	65	65	65	50
Min. dB HL	-10	-10	-10	-10	-10	-10	-10	-10	-10

6. APPLICABLE STANDARDS AND DIRECTIVES

- European Directive of Medical Devices 93/42/CEE
- ROHS Directive 2011/65/EU
- Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE).
- Regulation EC 1272/2008 on classification, labeling and packaging of substances and mixtures (REACH)
- Directiva 94/62/EC relativa a los envases y sus residuos, modificada por 2004/12/EC

- Quality System (EN ISO 13485:2012+AC:2012, EN ISO 9001:2008)
- Risk Management (EN ISO 14971:2012)
- Safety of medical equipment (EN 60601-1:2006+AC:2010+A11:2011+A1:2013+AC:2014)
- Electromagnetic Compatibility: (EN 60601-1-2:2015)
- Audiometría: EN 60645-1:2015
- Medical Device software: (EN 62304:2006+AC:2008+A1:2015)
- Usability: (EN 60601-1-6:2010+A1:2015 and EN 62366:2008+A1:2015)
- Biocompatibility: Biological evaluation of medical devices (EN ISO 10993-1:2009+AC:2010)
- Documentation and information: (EN 1041:2008, EN ISO 15223-1:2016)
- Ethernet: IEEE 802.3 10BaseT

OTHER STANDARDS

- EN ISO 389-1:2018 Air-Conduction Calibration (supra-aural headphones)
- EN ISO 389-3:2016 Bone-Conduction Calibration
- EN ISO 389-4:1998 Narrow-band Masking Noise Calibration
- EN ISO 8253-1:2010 Fundamental pure-tone air and bone conduction threshold audiometry

7. SYMBOLS

MANUFACTURER (manufacturer date, name and address of manufacturer)

PRODUCT'S REFERENCE

SERIAL NUMBER

LOT NUMBER

CE MARKING

STAND-BY

USB

ETHERNET

APPLICABLE PART BF

WARNING, RISK IDENTIFIED

WARNING, IDENTIFIED RISK

CONSULT INSTRUCTIONS FOR USE

DIRECT CURRENT

WASTE DISPOSAL OF ELECTRICAL / ELECTRONIC WASTE AGREEMENT TO THE WEEE DIRECTIVE

TEMPERATURE LIMIT

ATMOSPHERIC PRESSURE LIMITATION

HUMIDITY LIMITATION

KEEP DRY

FRAGILE, HANDLE WITH CARE

LIMITATION OF STACKING (MAXIMUM NUMBER OF BOXES)

THIS WAY UP PACKAGING

ANNEX 1: ELECTROMAGNETIC COMPATIBILITY

Guidance and manufacturer's declaration – electromagnetic emissions		
SIBELSOUND DUO is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment..		
Emissions test	Compliance	Electromagnetic environment - Guidance
RF (Radiated) emissions CISPR 11 (EN 55011)	Grupo 1 Clase B	
RF (Conducted) emissions CISPR 11 (EN 55011)	Grupo 1 Clase B	Not applicable when operating with batteries.
Harmonic emissions EN-IEC 61000-3-2	Class A	Not applicable when operating with batteries.
Voltage fluctuations / Flicker emissions EN-IEC 61000-3-3	Complies	Not applicable when operating with batteries.

Guidance and manufacturer's declaration – electromagnetic immunity			
SIBELSOUND DUO is intended for use in the electromagnetic environment specified below. The customer or the user of SIBELSOUND DUO should assure that it is used in such an environment. Below mentioned tests are applicable to enclosure. Additionally, test EN-IEC 61000-4-2 is applicable to input/output signal and patient connections.			
Immunity test	EN-IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN-IEC 61000-4-2	±8 kV contact ±2, ±4, ±8, ±15 kV, air	±8 kV contact ±2, ±4, ±8, ±15 kV, air	
Radiated RF EN-IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	3 V/m	
Proximity fields from RF communications EN-IEC 61000-4-3	(View following table).		Portable and mobile RF communications equipment should be used no closer than 30cm to any part of SIBELSOUND DUO, including cables. The distance is calculated according to: $d = \left[\frac{6}{E} \right] \sqrt{P}$ where P is the output power rating of the transmitter in watts (W), d is the recommended separation in meters (m) and E is the immunity level in V/m.
Power frequency (50 / 60 Hz) magnetic field EN-IEC 61000-4-8	30 A/m	30 A/m	Only applicable to devices sensitive to magnetic fields.
NOTE Ut is the a.c. mains voltage prior to application of the test level.			

Proximity fields EN-IEC 61000-4-3					
Frequency (MHz)	Band (MHz)	Modulation	Power (W)	Distance (m)	E (V/m)
385	380-390	Pulse 18Hz	1.8	0.3	27
450	430-470	FM D: ± 5 kHz 1 kHz sinus	2	0.3	28
710 745 780	704-787	Pulse 217Hz	0.2	0.3	9
810 870 930	800-960	Pulse 18Hz	2	0.3	28
1720 1845 1970	1700-1990	Pulse 217Hz	2	0.3	28
2450	2400-2570	Pulse 217Hz	2	0.3	28
5240 5500 5785	5100-5800	Pulse 217Hz	0.2	0.3	9

ANNEX 2: COMPLIANCE WITH DATA PROTECTION REGULATIONS

The aim of this section is to help users comply with applicable laws related to personal data processing, regulated by the **LOPD (Spanish Organic Law for Data Protection 15/1999, December 13th)**, and the new General Data Protection Regulation (**EU 2016/679 FROM THE EUROPEAN PARLIAMENT AND THE COUNCIL, April 27 2016**), on the protection of natural persons with regard to the processing of personal data and on the free movement of such data. The main points of such regulations are briefly described, specifying what needs to be taken into account in order to meet their legal requirements.

IMPORTANT NOTICE

- **According to the applicable laws, the user is the only person responsible for storing and processing their patients' data according to the law.**
- **Observing the recommendations included herein does not guarantee, in any case, the user's full compliance with data protection regulations.**

Introduction to Applicable Regulations

The applicable data protection regulation in Spain is, mainly, the Organic Law 15/1999 from December 13th, whereby public or private companies which hold automated files with private data, are required to implement a series of measures to guarantee the protection of such data, and to ensure that the persons concerned can fully exercise their rights.

Apart from the aforementioned organic law, there is another applicable regulation, ROYAL DECREE 1720/2007 from December 11th, which defines the implementation of technical and organizational measures necessary to guarantee the safety of automated files, data processing centers, premises, equipment, systems, programs, and the people involved in the automated processing of personal data.

Three security levels are established for personal data processing: basic, intermediate and high. Medical equipment stores data related to patients' health; therefore, they are classified as high level. Companies and professionals which handle personal data have, amongst others, the following obligations:

- Notification of the data files to the Spanish Data Protection Agency
- Adequacy of the data collection and processing to the LOPD (Spanish Organic Law on Data Protection)
- Drafting the Security Document

CONFIGURING THE SIBEL SOUND DUO PROTECTION

The audiometer includes the option to protect the unit with an access code or PIN. This option, which can be configured by the user, aims at protecting the equipment, and specifically the private data contained in it, against unauthorized access.

In order to comply with applicable regulations, the user must enable this option and configure their access code; it is the user's responsibility to provide this code to authorized persons (refer to section 3.5 in this manual).

DATA PROCESSING

This audiometer can transfer files with patients' data to a computer, which will then be used with the W50 Audiometry Software. This software also complies with the Spanish Data Protection Organic Law, as specified in its User Manual.

ANNEX 3: CALCULATION OF DIAGNOSTICS

1. 'MINISTRY' DIAGNOSIS

The percentage of hearing loss is determined from the tonal audiogram for the 500, 1000, 2000 and 3000 Hz frequencies, and according to the evaluation criteria of disabilities caused by hearing impairment, as specified by the Ministry of Labor and Social Affairs, BOE document no. 22 from January 26 2000.

A. Monoaural Hearing Loss

- Add up the hearing thresholds from the 4 frequencies for each ear.
- Hearing Loss (%) = Value from Table 1 corresponding to the sum of the thresholds. If the sum result is below 100 dB, it will not be classed as hearing loss.
- Average dB = (total thresholds) / 4.

A. Bilateral Hearing Loss

- Bilateral Loss (%) = $(5 * (\text{loss \% in the best ear}) + (\text{loss \% in the worst ear})) / 6$.

C. Diagnostic Classification of the Hearing Impairment

The MINISTRY diagnosis does not establish any diagnostic classification for hearing loss, although health specialist recommend implementing the following in the unit:

- NORMAL 0 to 10%
- SLIGHT HYPOACUSIA 11 to 30%
- MODERATE HYPOACUSIA 31 to 60%
- SEVERE HYPOACUSIA 61 to 100%

Table 1: Conversion of the estimated hearing level into a percentage of monoaural hearing loss

DSHL, dB	loss %	DSHL, dB	loss %	DSHL, dB	loss %	DSHL, dB	loss %
100	0.0	170	23.6	240	52.5	310	78.8
105	1.9	175	28.1	245	54.5	315	80.6
110	3.8	180	30.0	250	56.2	320	82.5
115	5.6	185	31.9	255	58.1	325	84.4
120	7.5	190	33.8	260	60.0	330	86.2
125	9.4	195	35.6	265	61.9	335	88.1
130	11.2	200	37.5	270	63.8	340	90.0
135	13.1	205	39.4	275	65.6	345	90.9
140	15.0	210	41.2	280	67.5	350	93.8
145	16.9	215	43.1	285	69.3	355	95.6
150	18.8	220	45.0	290	71.2	360	97.5
155	20.6	225	46.9	295	73.1	365	99.4
160	22.5	230	48.9	300	75.0	368	100.0
165	24.4	235	50.5	305	76.9		

2. 'COUNCIL' DIAGNOSIS

The percentage of hearing loss is determined from the tonal audiogram, and according to the 2nd formula of the Council of Physical Therapy. A percentage of hearing loss is applied to each ear at the 500, 1000, 2000 y 4000 Hz frequencies, according to the corresponding threshold.

A. Monoaural Hearing Loss

- In Table 2, look for the correspondence between the hearing thresholds for each frequency and their % of hearing loss.
- Hearing Loss (%) = Sum of the hearing losses for all four frequencies.
- Average dB = (total thresholds) / 4.

A. Bilateral Hearing Loss

- Bilateral Loss (%) = $(7 * (\text{loss \% in the best ear}) + (\text{loss \% in the worst ear})) / 8$.

C. Diagnostic Classification of the Hearing Impairment

The 'COUNCIL' diagnoses do not establish a hearing loss diagnostic classification. Therefore, the 'MINISTRY' diagnosis classification is used.

Table 2: Council of Physical Therapy, 2nd Formula

Loss in dB	500 (Hz)	1000 (Hz)	2000 (Hz)	4000 (Hz)
5				
10	0.2	0.3	0.4	0.1
15	0.5	0.9	1.3	0.3
20	1.1	2.1	2.9	0.9
25	1.8	3.6	4.9	1.7
30	2.6	5.4	7.3	2.7
35	3.7	7.7	9.8	3.8
40	4.9	10.2	12.9	5.0
45	6.3	13.0	17.3	6.4
50	7.9	15.7	22.4	8.0
55	9.6	19.0	25.7	9.7
60	11.3	21.5	28.0	11.2
65	12.8	23.5	30.2	12.5
70	13.8	25.5	32.2	13.5
75	14.6	27.2	34.4	14.2
80	14.8	28.8	35.8	14.6
85	14.9	29.8	37.5	14.8
90	15.0	29.9	39.2	14.9
95	15.0	30.0	40.0	15.0
100				

3. DIAGNOSIS IN ACCORDANCE WITH THE MEXICAN SOCIAL SECURITY INSTITUTE

The percentage of hearing loss is determined from the tonal audiogram for the 500, 1000, 2000 and 3000 Hz frequencies, and according to the Mexican Social Security Institute.

A. Monoaural Hearing Loss

- Hearing Loss (%) = (total thresholds / 4) * 0.8.
- Average dB = (total thresholds / 4).

A. Bilateral Hearing Loss

- Bilateral Loss (%) = (7 * (loss % in the best ear) + (loss % in the worst ear)) / 8.

C. Diagnostic Classification of the Hearing Impairment

The 'COUNCIL' diagnoses do not establish a hearing loss diagnostic classification. Therefore, the 'MINISTRY' diagnosis classification is used.

4. ELI INDEX

The ELI Index (Early Loss Index) is calculated from the 4000 Hz frequency, following the Specific Health Surveillance Protocol (Spanish Interterritorial Council of the National Health Service). The percentages of hearing loss are not indicated.

A. ELI Index

ELI Index = Threshold at 4000 Hz – Table 3 value according to age and gender

B. Diagnostic Classification of the Hearing Impairment

In Table 4, acoustic traumas are classified with an increasing A-B-C-D-E scale, from greater to lesser hearing capacity.

Table 3: Table for calculating the ELI Index, Presbycusis correction, at 4,000 Hz, dB

Age	Women	Men
25	0	0
30	2	3
35	3	7
40	5	11
45	8	15
50	12	20
55	15	26
60	17	32
65	18	38

Table 4: Table for classifying acoustic traumas

Cóorrected Audiometric Loss, dB	Degree of ELI	Classification
<8	A	Normal excellent
8-14	B	Normal good
15-22	C	Normal
23-29	D	Suspicion of deafness
>30	E	Clear indication of deafness

5. SAL INDEX

The SAL Index (Speech Average Loss) assesses conversational frequencies (500, 1000 and 2000 Hz), and is defined as the arithmetic mean of hearing loss, in decibels, at those frequencies

A. Monoaural and Bilateral Hearing Loss

They are calculated in the same way as the MINISTRY diagnosis.

B. Diagnostic Classification of the Hearing Impairment

As per table 5.

When the best ear presents a decline of <16 dB and the worst ear >30 dB, it is considered that there are alterations which are excessively asymmetric, not considered in table 5. In this case, the following text appears:

'The deafness detected is excessively asymmetric to be evaluated using the SAL table'

Table 5: Evaluation and meaning of the SAL Index

Degree of Sal	dB	Class name	Characteristics
A	16 worse ear	Normal	Both ears are within normal limits, with no difficulty in quiet conversations.
B	16-30 one of the ears	Close to normal	Has difficulty with quiet conversations but nothing more.
C	31-45 better ear	Slight Auditory Loss	Has difficulties in normal conversations, but not if people speak up.
D	46-60 better ear	Serious Auditory Loss	Has difficulties even when people speak up.
E	61-90 better ear	Serious Auditory Loss	Can only hear a conversation using a hearing aid.
F	91-99 better ear	Profound Auditory Loss	Cannot hear a conversation even with a hearing aid.
G	> 100 better ear	Totally deaf in both ears	Can hear no sound at all.

6. KLOCKHOOF DIAGNOSIS (Modified by Lavoro Clinic in Milan)

This diagnosis evaluates the 500, 1000, 2000 and 4000 Hz frequencies. The following criteria apply in accordance with Table 6:

Table 6: KLOCKHOFF Classification

NORMAL	The threshold does not exceed 25 dB at any frequency.		
ISOLATED ACOUSTIC TRAUMA	No conversational loss	LIGHT	The threshold at the 4000 frequency is below 55 dB and below 25 dB in the other frequencies.
		ADVANCED	The threshold at the 4000 frequency is above 55 dB and below 25 dB at the other frequencies.
GLOBAL ACOUSTIC TRAUMA	Conversational loss	LIGHT	One or more frequencies are not affected (are below 25 dB).
		MODERATE	All frequencies are affected (above 25 dB) but none is above 55 dB.
		ADVANCED	All frequencies are affected (above 25 dB) and one above 55 dB.
OTHER LOSS not due to exposure to sound.			

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All rights reserved.